

# PRODUCT CATALOG



***Implantech***<sup>®</sup>  
Superior Patient Aesthetics

# Superior Facial & Body Aesthetics

Designed by prominent surgeons, Implantech's facial implants provide the youthful, cosmetic and reconstructive solutions desired by your patients. Count on our comprehensive line of silicone facial implants for industry-leading innovation and effective, pleasing outcomes.

If no stock product is appropriate, tap into our patient-specific resources. For any part of the face or body, we faithfully reproduce silicone implants from CT scans with our 3D Accuscan® service. Or, carve your own from silicone blocks of different sizes, shapes and durometers.

For scar management, our Cimeosil® products are specifically formulated for the topical management of keloids and hypertrophic scars. Our Cimeosil® Scar and Laser Gel also helps with any scarring or redness associated with laser surgery, dermabrasion and chemical peels. Made of stretchable, fabric-covered silicone gel, our Gelzone® Wraps combine compression healing, musculoskeletal support and scar management.

You can also depend on Implantech for sheeting. We offer silicone sheeting for short term or long term applications.

Get the service, quality, innovation and choices you and your patients deserve. Welcome to Implantech!

# Table of Contents

|                                                                |    |
|----------------------------------------------------------------|----|
| <b><i>Silicone Facial Implants</i></b> .....                   | 2  |
| Chin .....                                                     | 3  |
| Chin Sizing Guide .....                                        | 7  |
| Mandibular Angle.....                                          | 9  |
| Submalar .....                                                 | 11 |
| Midfacial.....                                                 | 12 |
| Tear Trough.....                                               | 13 |
| Temporal .....                                                 | 14 |
| Nasal .....                                                    | 15 |
| <br><b><i>Silicone Carving Blocks</i></b> .....                | 17 |
| Rectangular Carving Blocks .....                               | 18 |
| Contoured Carving Blocks.....                                  | 19 |
| <br><b><i>3D Accuscan Patient-Specific® Implants</i></b> ..... | 21 |
| Patient-Specific Resources (3D Accuscan®) .....                | 22 |
| <br><b><i>Silicone Sheeting</i></b> .....                      | 23 |
| AlliedSil™ Silicone Sheeting (Non-Reinforced).....             | 24 |
| <br><b><i>Scar Management</i></b> .....                        | 25 |
| Cimeosil® Scar and Laser Gel .....                             | 26 |
| Cimeosil® Gel Sheeting.....                                    | 26 |
| <br><b><i>Gelzone® Wrap</i></b> .....                          | 28 |
| Gelzone® Wrap.....                                             | 28 |
| <br><b><i>Ordering Information</i></b> .....                   | 29 |

# SILICONE FACIAL IMPLANTS



## ***All Supplied Sterile***

Crafted with the highest quality of implant grade silicone, our implants are ideal for augmenting chin, midfacial, mandibular, nasal and other facial areas.

Choose from over 35 styles and variations of facial implants for a long term and symmetric look, including:


- Our core line of trusted, time-tested silicone implants
- Conform™ silicone implants with a grid-like backing that readily adapts to bony contours

## Conform™ Extended Anatomical Chin

Evolution of all of the technology and design of our previous chin implants. Incorporates the Conform™ grid pattern for flexible placement and easier contouring along the underlying bone.

| Cat. No.    | Size    | A            | B    | C    | D    | E   | F    |
|-------------|---------|--------------|------|------|------|-----|------|
| CEAC-S      | Small   | 5.0          | 1.00 | 0.43 | 8.7  | 3.5 | 0.22 |
| CEAC-M      | Medium  | 5.6          | 1.20 | 0.50 | 10.5 | 4.2 | 0.25 |
| CEAC-L      | Large   | 6.2          | 1.35 | 0.60 | 11.5 | 4.5 | 0.31 |
| CEAC-XL     | X-Large | 6.8          | 1.50 | 0.75 | 12.5 | 5.0 | 0.33 |
| CEAC-SIZ-NS | Sizers  | Complete Set |      |      |      |     |      |

*Dimensions in centimeters*




## Extended Anatomical Chin

Full lateral arms blend into mandible. Can be used with Glasgold Wafer™ (GLE-2) for additional anterior augmentation.

| Cat. No.   | Size    | A            | B   | C    | D    | E   | F   |
|------------|---------|--------------|-----|------|------|-----|-----|
| EAC-S      | Small   | 5.5          | 1.0 | 0.45 | 9.0  | 3.7 | 0.2 |
| EAC-M      | Medium  | 5.6          | 1.2 | 0.50 | 9.5  | 4.2 | 0.4 |
| EAC-L      | Large   | 6.8          | 1.4 | 0.60 | 11.0 | 4.8 | 0.4 |
| EAC-XL     | X-Large | 6.0          | 1.4 | 0.75 | 11.0 | 5.1 | 0.5 |
| EAC-SIZ-NS | Sizers  | Complete Set |     |      |      |     |     |

*Dimensions in centimeters*




## Terino Extended Anatomical™ Chin


Lateral arms are thinner and broader than Extended Anatomical Chin for less palpable edges. Profile is slightly squarer for a stronger jawline. Can be used with Glasgold Wafer™ (GLE-2) for additional anterior augmentation.

| Cat. No.    | Size    | A            | B   | C    | D    | E   | F   |
|-------------|---------|--------------|-----|------|------|-----|-----|
| TEAC-S      | Small   | 5.7          | 1.0 | 0.38 | 11.0 | 4.8 | 0.3 |
| TEAC-M      | Medium  | 5.7          | 1.2 | 0.46 | 11.3 | 5.0 | 0.3 |
| TEAC-L      | Large   | 5.7          | 1.3 | 0.58 | 11.5 | 5.2 | 0.3 |
| TEAC-XL     | X-Large | 5.7          | 1.4 | 0.64 | 11.5 | 5.2 | 0.4 |
| TEAC-SIZ-NS | Sizers  | Complete Set |     |      |      |     |     |

*Dimensions in centimeters*



## CHIN IMPLANTS




### Flowers Mandibular Glove®

Posterior implant surface matches the anatomical variations of the anterior surface of the mandible. Provides anterior chin augmentation while slightly tilting the chin mound upward. Blends in laterally, with no appreciable lateral augmentation.

| Cat. No.   | Size    | A            | B    | C    | D    | E    | F   |
|------------|---------|--------------|------|------|------|------|-----|
| FMG-VXS    | X-Small | 6.00         | 1.75 | 0.50 | 9.0  | 3.65 | 0.2 |
| FMG-VS     | Small   | 6.02         | 1.75 | 0.60 | 9.0  | 3.75 | 0.3 |
| FMG-VM     | Medium  | 6.06         | 1.75 | 0.80 | 10.0 | 4.41 | 0.3 |
| FMG-VL     | Large   | 6.10         | 1.75 | 1.00 | 10.0 | 4.61 | 0.4 |
| FMG-VXL    | X-Large | 6.30         | 1.75 | 1.40 | 10.8 | 5.40 | 0.5 |
| FMG-SIZ-NS | Sizers  | Complete Set |      |      |      |      |     |

*Dimensions in centimeters*




### Vertical Lengthening Chin

**Designed by Barry L. Eppler, MD, DMD**

Designed to address a three-dimensional chin deficiency. Implant secures to the lower edge of the anterior jawline to provide both a downward 45-degree augmentation and a horizontal projection. Creates a more elongated face and triangular jawline-chin shape.

| Cat. No.   | Size   | A            | B   | C    | D    | E   | F    | G    |
|------------|--------|--------------|-----|------|------|-----|------|------|
| VLC-S      | Small  | 6.5          | 1.4 | 0.79 | 9.0  | 3.6 | 0.52 | 0.55 |
| VLC-M      | Medium | 7.3          | 1.6 | 1.00 | 10.5 | 4.7 | 0.50 | 0.60 |
| VLC-L      | Large  | 7.8          | 1.8 | 1.20 | 12.0 | 5.2 | 0.60 | 0.70 |
| VLC-SIZ-NS | Sizers | Complete Set |     |      |      |     |      |      |

*Dimensions in centimeters*



### Mandibular Pre Jowl®

Provides augmentation in pre-jowl sulcus with no additional anterior chin augmentation.

| Cat. No.   | Size    | A            | B   | C    | D    | E   | F   |
|------------|---------|--------------|-----|------|------|-----|-----|
| MPJ-S      | Small   | 6.0          | 1.0 | 0.10 | 10.9 | 4.2 | 0.3 |
| MPJ-M      | Medium  | 6.0          | 1.0 | 0.10 | 11.0 | 4.3 | 0.4 |
| MPJ-L      | Large   | 6.0          | 1.0 | 0.10 | 11.2 | 4.3 | 0.5 |
| MPJ-XL     | X-Large | 6.0          | 1.0 | 0.10 | 12.5 | 4.6 | 0.6 |
| MPJ-SIZ-NS | Sizers  | Complete Set |     |      |      |     |     |


*Dimensions in centimeters*

## Mandibular Pre Jowl Chin™

Provides anterior chin augmentation with a subtle increase in lateral fullness precisely positioned to augment the pre-jowl sulcus. Posterior implant surface allows for a pronounced symphyseal ridge. Can be used with Glasgold Wafer™ (GLM-2) for additional anterior augmentation.

| Cat. No.    | Size    | A            | B   | C    | D    | E   | F    |
|-------------|---------|--------------|-----|------|------|-----|------|
| MPJC-S      | Small   | 6.0          | 1.5 | 0.65 | 10.5 | 4.0 | 0.33 |
| MPJC-M      | Medium  | 6.0          | 1.5 | 0.80 | 10.8 | 4.3 | 0.40 |
| MPJC-L      | Large   | 6.0          | 1.5 | 0.95 | 11.0 | 4.8 | 0.50 |
| MPJC-XL     | X-Large | 6.0          | 1.5 | 1.10 | 11.0 | 5.0 | 0.56 |
| MPJC-SIZ-NS | Sizers  | Complete Set |     |      |      |     |      |

*Dimensions in centimeters*




## Terino Square Chin - Style I

Provides anterior chin augmentation with a square profile.

| Cat. No.    | Size    | A            | B   | C    | D    | E   | F   | G   | H    |
|-------------|---------|--------------|-----|------|------|-----|-----|-----|------|
| TSCI-S      | Small   | 6.6          | 1.3 | 0.60 | 9.8  | 4.1 | 0.4 | 4.0 | 0.70 |
| TSCI-M      | Medium  | 6.5          | 1.4 | 0.70 | 10.9 | 4.8 | 0.3 | 4.2 | 0.80 |
| TSCI-L      | Large   | 6.6          | 1.6 | 0.80 | 12.0 | 5.4 | 0.3 | 4.6 | 0.93 |
| TSCI-XL     | X-Large | 6.9          | 1.6 | 1.10 | 12.8 | 5.8 | 0.4 | 4.2 | 1.20 |
| TSCI-SIZ-NS | Sizers  | Complete Set |     |      |      |     |     |     |      |

*Dimensions in centimeters*





## Terino Square Chin - Style II

A more strongly defined square profile than Style I.

| Cat. No.     | Size    | A            | B   | C    | D    | E   | F   | G   | H    |
|--------------|---------|--------------|-----|------|------|-----|-----|-----|------|
| TSCII-S      | Small   | 7.8          | 1.5 | 0.47 | 11.5 | 4.3 | 0.4 | 4.0 | 0.74 |
| TSCII-M      | Medium  | 8.0          | 1.5 | 0.50 | 12.8 | 5.0 | 0.4 | 5.4 | 0.80 |
| TSCII-L      | Large   | 9.4          | 1.5 | 0.55 | 13.3 | 4.9 | 0.4 | 5.1 | 0.91 |
| TSCII-XL     | X-Large | 9.0          | 1.5 | 0.65 | 15.0 | 6.2 | 0.4 | 5.3 | 1.00 |
| TSCII-SIZ-NS | Sizers  | Complete Set |     |      |      |     |     |     |      |

*Dimensions in centimeters*






### Anatomical Chin

Central chin augmentation. Short lateral extensions blend into mandible.

| Cat. No.  | Size    | A            | B   | C    | D   | E   |
|-----------|---------|--------------|-----|------|-----|-----|
| AC-S      | Small   | 4.7          | 1.2 | 0.60 | 6.0 | 2.1 |
| AC-M      | Medium  | 5.3          | 1.2 | 0.90 | 6.6 | 2.5 |
| AC-L      | Large   | 5.1          | 1.2 | 1.00 | 6.2 | 2.5 |
| AC-XL     | X-Large | 5.3          | 1.2 | 1.20 | 6.6 | 2.8 |
| AC-SIZ-NS | Sizers  | Complete Set |     |      |     |     |

*Dimensions in centimeters*




### Curvilinear Silicone Chin

Central chin augmentation with no lateral component.

| Cat. No.   | Size    | A            | B   | C    | D   |
|------------|---------|--------------|-----|------|-----|
| CSC-XS     | X-Small | 3.1          | 1.1 | 0.45 | 3.0 |
| CSC-S      | Small   | 3.2          | 1.4 | 0.65 | 3.3 |
| CSC-M      | Medium  | 3.6          | 1.4 | 0.65 | 4.0 |
| CSC-L      | Large   | 4.4          | 1.4 | 0.65 | 4.8 |
| CSC-XL     | X-Large | 4.0          | 1.5 | 0.70 | 4.5 |
| CSC-SIZ-NS | Sizers  | Complete Set |     |      |     |

*Dimensions in centimeters*



### Glasgold Wafer™

Provides 2mm additional anterior projection for Extended Anatomical Chin, Terino Extended Anatomical™ Chin (GLE-2) and Mandibular Pre Jowl Chin™ (GLM-2).

| Cat. No. | Size | A   | B   | C   |
|----------|------|-----|-----|-----|
| GLE-2    | 2mm  | 3.6 | 1.1 | 0.2 |
| GLM-2    | 2mm  | 3.6 | 1.1 | 0.2 |

*Dimensions in centimeters*

### Tips on Choosing the Appropriate Implantech Chin Implant

Table 1. DESIRED OUTCOME – Compare the three dimensions of projection, by style and size

| Anterior Projection                                 | Cat. No.    | X-Small                                           | Small | Medium | Large | X-Large |
|-----------------------------------------------------|-------------|---------------------------------------------------|-------|--------|-------|---------|
| Terino Extended Anatomical™ Chin                    | TEAC        |                                                   | 0.38  | 0.46   | 0.58  | 0.64    |
| Conform™ Extended Anatomical Chin                   | CEAC        |                                                   | 0.43  | 0.50   | 0.60  | 0.75    |
| * Extended Anatomical Chin                          | EAC         |                                                   | 0.45  | 0.50   | 0.60  | 0.75    |
| * Curvilinear Silicone Chin (no lateral projection) | CSC         | 0.45                                              | 0.65  | 0.65   | 0.65  | 0.70    |
| Flowers Mandibular Glove®                           | FMG         | 0.50                                              | 0.60  | 0.80   | 1.00  | 1.40    |
| * Anatomical Chin (no lateral projection)           | AC          |                                                   | 0.60  | 0.90   | 1.00  | 1.20    |
| Mandibular Pre Jowl Chin™                           | MPJC        |                                                   | 0.65  | 0.80   | 0.95  | 1.10    |
| Vertical Lengthening Chin                           | VLC         |                                                   | 0.79  | 1.00   | 1.20  |         |
| Glasgold Wafer™                                     | GLE-2/GLM-2 | Provides additional 0.20 cm to EAC, TEAC and MPJC |       |        |       |         |
| Lateral Projection                                  | Cat. No.    | X-Small                                           | Small | Medium | Large | X-Large |
| Conform™ Extended Anatomical Chin                   | CEAC        |                                                   | 0.22  | 0.25   | 0.31  | 0.33    |
| Flowers Mandibular Glove®                           | FMG         | 0.20                                              | 0.30  | 0.30   | 0.40  | 0.50    |
| * Extended Anatomical Chin                          | EAC         |                                                   | 0.20  | 0.40   | 0.40  | 0.50    |
| Terino Extended Anatomical™ Chin                    | TEAC        |                                                   | 0.30  | 0.30   | 0.30  | 0.40    |
| Mandibular Pre Jowl® (no ant. proj. or vert. ht.)   | MPJ         |                                                   | 0.30  | 0.40   | 0.50  | 0.60    |
| Mandibular Pre Jowl Chin™                           | MPJC        |                                                   | 0.33  | 0.40   | 0.50  | 0.56    |
| Vertical Lengthening Chin                           | VLC         |                                                   | 0.52  | 0.50   | 0.60  |         |
| Vertical Height                                     | Cat. No.    | X-Small                                           | Small | Medium | Large | X-Large |
| Terino Extended Anatomical™ Chin                    | TEAC        |                                                   | 1.00  | 1.20   | 1.30  | 1.40    |
| * Extended Anatomical Chin                          | EAC         |                                                   | 1.00  | 1.20   | 1.40  | 1.40    |
| Conform™ Extended Anatomical Chin                   | CEAC        |                                                   | 1.00  | 1.20   | 1.35  | 1.50    |
| * Anatomical Chin                                   | AC          |                                                   | 1.20  | 1.20   | 1.20  | 1.20    |
| * Curvilinear Silicone Chin                         | CSC         | 1.10                                              | 1.40  | 1.40   | 1.40  | 1.50    |
| Vertical Lengthening Chin                           | VLC         |                                                   | 1.40  | 1.60   | 1.80  |         |
| Mandibular Pre Jowl Chin™                           | MPJC        |                                                   | 1.50  | 1.50   | 1.50  | 1.50    |
| Flowers Mandibular Glove®                           | FMG         | 1.75                                              | 1.75  | 1.75   | 1.75  | 1.75    |

Dimensions in centimeters

## Other Desired Characteristics

### Square Appearance

Terino Extended Anatomical™ Chin (TEAC)  
Terino Square Chin-Style I (TSCI)  
Terino Square Chin-Style II (TSCII)

### Slight tilt upwards

Flowers Mandibular Glove® (FMG)

### Vertical Dimension

Vertical Lengthening Chin (VLC)

### Stability

Implants with holes for integration/fixation

Curvilinear Silicone Chin (CSC)

Mandibular Pre Jowl® (MPJ)

Mandibular Pre Jowl Chin™ (MPJC)

Conforming (more flexible)

Conform™ Extended Anatomical Chin (CEAC)

## Tips on Choosing the Appropriate Implantech Chin Implant

Table 2. ANTERIOR PROJECTION – Regardless of style, from least to greatest

| Implant Style                     | Cat. No. | Size    | Anterior Projection |
|-----------------------------------|----------|---------|---------------------|
| Terino Extended Anatomical™ Chin  | TEAC-S   | Small   | 0.38                |
| Conform™ Extended Anatomical Chin | CEAC-S   | Small   | 0.43                |
| * Extended Anatomical Chin        | EAC-S    | Small   | 0.45                |
| * Curvilinear Silicone Chin       | CSC-XS   | X-Small | 0.45                |
| Terino Extended Anatomical™ Chin  | TEAC-M   | Medium  | 0.46                |
| Conform™ Extended Anatomical Chin | CEAC-M   | Medium  | 0.50                |
| * Extended Anatomical Chin        | EAC-M    | Medium  | 0.50                |
| Flowers Mandibular Glove®         | FMG-VXS  | X-Small | 0.50                |
| Terino Extended Anatomical™ Chin  | TEAC-L   | Large   | 0.58                |
| Conform™ Extended Anatomical Chin | CEAC-L   | Large   | 0.60                |
| * Extended Anatomical Chin        | EAC-L    | Large   | 0.60                |
| Flowers Mandibular Glove®         | FMG-VS   | Small   | 0.60                |
| * Anatomical Chin                 | AC-S     | Small   | 0.60                |
| Terino Extended Anatomical™ Chin  | TEAC-XL  | X-Large | 0.64                |
| Mandibular Pre Jowl Chin™         | MPJC-S   | Small   | 0.65                |
| * Curvilinear Silicone Chin       | CSC-S    | Small   | 0.65                |
| * Curvilinear Silicone Chin       | CSC-M    | Medium  | 0.65                |
| * Curvilinear Silicone Chin       | CSC-L    | Large   | 0.65                |
| * Curvilinear Silicone Chin       | CSC-XL   | X-Large | 0.70                |
| Conform™ Extended Anatomical Chin | CEAC-XL  | X-Large | 0.75                |
| * Extended Anatomical Chin        | EAC-XL   | X-Large | 0.75                |
| Vertical Lengthening Chin         | VLC-S    | Small   | 0.79                |
| Flowers Mandibular Glove®         | FMG-VM   | Medium  | 0.80                |
| Mandibular Pre Jowl Chin™         | MPJC-M   | Medium  | 0.80                |
| * Anatomical Chin                 | AC-M     | Medium  | 0.90                |
| Mandibular Pre Jowl Chin™         | MPJC-L   | Large   | 0.95                |
| Vertical Lengthening Chin         | VLC-M    | Medium  | 1.00                |
| Flowers Mandibular Glove®         | FMG-VL   | Large   | 1.00                |
| * Anatomical Chin                 | AC-L     | Large   | 1.00                |
| Mandibular Pre Jowl Chin™         | MPJC-XL  | X-Large | 1.10                |
| Vertical Lengthening Chin         | VLC-L    | Large   | 1.20                |
| * Anatomical Chin                 | AC-XL    | X-Large | 1.20                |
| Flowers Mandibular Glove®         | FMG-VXL  | X-Large | 1.40                |

Dimensions in centimeters


Chin Sizing Guide suggested by Steven J. Pearlman, MD, FACS

## Conform™ Mandibular Angle

Compared to the Lateral Mandibular Angle (TMA), this implant's shortened ramus, longer mandibular dimensions and slightly more lateral projection put more emphasis on the angle prominence and mandibular definition. Conform™ grid pattern readily adapts to the gonial angle. Specify Catalog No. for "Left (L)" or "Right (R)". Single-packed (not in pairs). Save by purchasing a pair.

| Cat. No.   | A      | B            | C    | D   | E   | F   | G    |
|------------|--------|--------------|------|-----|-----|-----|------|
| CMA-8L     | 3.70   | 0.40         | 4.70 | 3.0 | 3.0 | 0.8 | 0.40 |
| CMA-8R     | 3.70   | 0.40         | 4.70 | 3.0 | 3.0 | 0.8 | 0.40 |
| CMA-10L    | 3.70   | 0.40         | 4.80 | 3.0 | 3.0 | 1.0 | 0.40 |
| CMA-10R    | 3.70   | 0.40         | 4.80 | 3.0 | 3.0 | 1.0 | 0.40 |
| CMA-12L    | 3.70   | 0.40         | 4.90 | 3.0 | 3.0 | 1.2 | 0.40 |
| CMA-12R    | 3.70   | 0.40         | 4.90 | 3.0 | 3.0 | 1.2 | 0.40 |
| CMA-14L    | 3.70   | 0.40         | 5.00 | 3.0 | 3.0 | 1.4 | 0.40 |
| CMA-14R    | 3.70   | 0.40         | 5.00 | 3.0 | 3.0 | 1.4 | 0.40 |
| CMA-SIZ-NS | Sizers | Complete Set |      |     |     |     |      |

Dimensions in centimeters




## Widening Mandibular Angle

Designed by Barry L. Eppley, MD, DMD

Features less lateral projection than the Conform™ Mandibular Angle (CMA) and a more defined angle prominence. Extended mandible and shortened ramus edges help enable a more ideal fit to jaw angle. Specify Catalog No. for "Left (L)" or "Right (R)". Single-packed (not in pairs). Save by purchasing a pair.

| Cat. No.   | A      | B            | C   |
|------------|--------|--------------|-----|
| WMA-5L     | 3.2    | 0.5          | 4.0 |
| WMA-5R     | 3.2    | 0.5          | 4.0 |
| WMA-7L     | 3.3    | 0.7          | 4.4 |
| WMA-7R     | 3.3    | 0.7          | 4.4 |
| WMA-10L    | 3.6    | 1.0          | 4.6 |
| WMA-10R    | 3.6    | 1.0          | 4.6 |
| WMA-SIZ-NS | Sizers | Complete Set |     |

Dimensions in centimeters




## Vertical Mandibular Angle

Designed by Barry L. Eppley, MD, DMD


An innovative jaw implant design that creates a distinct jaw angle shape below the patient's natural jawline. Vertically lengthens the short ramus to more clearly define the jaw angle. Specify Catalog No. for "Left (L)" or "Right (R)". Single-packed (not in pairs). Save by purchasing a pair.

| Cat. No.   | A      | B            | C   | D   |
|------------|--------|--------------|-----|-----|
| VMA-7L     | 2.5    | 0.4          | 4.5 | 0.7 |
| VMA-7R     | 2.5    | 0.4          | 4.5 | 0.7 |
| VMA-11L    | 3.0    | 0.5          | 5.0 | 1.1 |
| VMA-11R    | 3.0    | 0.5          | 5.0 | 1.1 |
| VMA-13L    | 3.6    | 0.5          | 5.5 | 1.3 |
| VMA-13R    | 3.6    | 0.5          | 5.5 | 1.3 |
| VMA-SIZ-NS | Sizers | Complete Set |     |     |

Dimensions in centimeters



## MANDIBULAR ANGLE IMPLANTS




### Lateral Mandibular Angle™

Lateral augmentation of mandibular angle (dimension F). Specify Catalog No. for "Left (L)" or "Right (R)." Lateral Mandibular Implants are single-packed (not in pairs).

| Cat. No.     | A      | B            | C   | D   | E   | F   | G   |
|--------------|--------|--------------|-----|-----|-----|-----|-----|
| TMA-4L/8     | 4.5    | 0.4          | 4.5 | 3.5 | 3.0 | 0.8 | 0.4 |
| TMA-4R/8     | 4.5    | 0.4          | 4.5 | 3.5 | 3.0 | 0.8 | 0.4 |
| TMA-4L/10    | 4.5    | 0.4          | 4.5 | 3.5 | 3.0 | 1.0 | 0.4 |
| TMA-4R/10    | 4.5    | 0.4          | 4.5 | 3.5 | 3.0 | 1.0 | 0.4 |
| TMA-4L/12    | 4.5    | 0.4          | 4.5 | 3.5 | 3.0 | 1.2 | 0.4 |
| TMA-4R/12    | 4.5    | 0.4          | 4.5 | 3.5 | 3.0 | 1.2 | 0.4 |
| TMA-L-SIZ-NS | Sizers | Complete Set |     |     |     |     |     |

*Dimensions in centimeters*



### Posterior Mandibular Angle™

Augments posterior mandibular angle (dimension B) with minor lateral augmentation (dimension F). Specify Catalog No. for "Left (L)" or "Right (R)." Posterior Mandibular Implants are single-packed (not in pairs).

| Cat. No. | A   | B   | C   | D   | E   | F    | G   |
|----------|-----|-----|-----|-----|-----|------|-----|
| TMA-4L   | 4.5 | 0.4 | 4.5 | 3.5 | 3.0 | 0.36 | 0.4 |
| TMA-4R   | 4.5 | 0.4 | 4.5 | 3.5 | 3.0 | 0.36 | 0.4 |
| TMA-6L   | 4.5 | 0.6 | 4.5 | 3.5 | 3.0 | 0.36 | 0.4 |
| TMA-6R   | 4.5 | 0.6 | 4.5 | 3.5 | 3.0 | 0.36 | 0.4 |
| TMA-8L   | 4.5 | 0.8 | 4.5 | 3.5 | 3.0 | 0.36 | 0.4 |
| TMA-8R   | 4.5 | 0.8 | 4.5 | 3.5 | 3.0 | 0.36 | 0.4 |
| TMA-10L  | 4.5 | 1.0 | 4.5 | 3.5 | 3.0 | 0.36 | 0.4 |
| TMA-10R  | 4.5 | 1.0 | 4.5 | 3.5 | 3.0 | 0.36 | 0.4 |


*Dimensions in centimeters*

## Conform Binder Submalar®

Making it more flexible than the standard silicone implant, the Conform™ grid pattern readily adapts to the underlying facial structure. Augments the submalar region for a more youthful appearance. Sold in pairs.

| Cat. No.    | Size    | A            | B   | C   |
|-------------|---------|--------------|-----|-----|
| CBSM-S      | Small   | 4.4          | 1.5 | 0.4 |
| CBSM-M      | Medium  | 4.9          | 1.8 | 0.5 |
| CBSM-L      | Large   | 5.4          | 1.9 | 0.6 |
| CBSM-XL     | X-Large | 5.9          | 2.1 | 0.7 |
| CBSM-SIZ-NS | Sizers  | Complete Set |     |     |

Dimensions in centimeters




## Binder Submalar®

Elevates and repositions the overlying midfacial soft tissue in an anterior and superior position for a more youthful appearance. Adds lower midfacial fullness while softening the nasolabial fold. Sold in pairs.

| Cat. No.   | Size    | A            | B    | C    |
|------------|---------|--------------|------|------|
| BSM-S      | Small   | 4.03         | 1.62 | 0.45 |
| BSM-M      | Medium  | 4.48         | 1.82 | 0.50 |
| BSM-L      | Large   | 4.97         | 2.02 | 0.56 |
| BSM-XL     | X-Large | 5.33         | 2.23 | 0.61 |
| BSM-SIZ-NS | Sizers  | Complete Set |      |      |

Dimensions in centimeters




## Binder Submalar® II


Same aesthetic benefit as Binder Submalar®, but designed to fit more inferiorly into the submalar space to provide natural looking volume in patients with high cheekbones and sunken cheeks. Sold in pairs.

| Cat. No.     | Size    | A            | B    | C    |
|--------------|---------|--------------|------|------|
| BSMII-S      | Small   | 3.95         | 2.14 | 0.40 |
| BSMII-M      | Medium  | 4.38         | 2.26 | 0.45 |
| BSMII-L      | Large   | 4.88         | 2.54 | 0.54 |
| BSMII-XL     | X-Large | 5.22         | 2.77 | 0.54 |
| BSMII-SIZ-NS | Sizers  | Complete Set |      |      |

Dimensions in centimeters



## MIDFACIAL IMPLANTS



### Conform Terino Malar Shell®


Augments the volume deficient midface for a more aesthetic, youthful and balanced look. The Conform™ design is flexible and enables the implant to more easily adapt to the underlying midfacial structure. Sold in pairs.

| Cat. No.    | Size    | A            | B   | C   |
|-------------|---------|--------------|-----|-----|
| CTMS-S4     | Small   | 4.8          | 2.3 | 0.4 |
| CTMS-M4     | Medium  | 5.1          | 2.8 | 0.4 |
| CTMS-L4     | Large   | 5.2          | 3.3 | 0.4 |
| CTMS-XL4    | X-Large | 5.7          | 3.4 | 0.4 |
| CTMS-SIZ-NS | Sizers  | Complete Set |     |     |

*Dimensions in centimeters*


### Terino Malar Shell®

Broad surface area with fine, non-palpable edges for a smooth transition between the implant and underlying bone. Ideal for augmenting patients with malar hypoplasia. Sold in pairs.



| Cat. No.   | Size    | A                              | B   | C   |
|------------|---------|--------------------------------|-----|-----|
| TMS-S3     | Small   | 4.8                            | 2.3 | 0.3 |
| TMS-S4     | Small   | 4.8                            | 2.3 | 0.4 |
| TMS-M3     | Medium  | 5.1                            | 2.8 | 0.3 |
| TMS-M4     | Medium  | 5.1                            | 2.8 | 0.4 |
| TMS-M5     | Medium  | 5.1                            | 2.8 | 0.5 |
| TMS-L3     | Large   | 5.2                            | 3.3 | 0.3 |
| TMS-L4     | Large   | 5.2                            | 3.3 | 0.4 |
| TMS-L5     | Large   | 5.2                            | 3.3 | 0.5 |
| TMS-L6     | Large   | 5.2                            | 3.3 | 0.6 |
| TMS-XL4    | X-Large | 6.0                            | 3.7 | 0.4 |
| TMS-XL5    | X-Large | 6.0                            | 3.7 | 0.5 |
| TMS-XL6    | X-Large | 6.0                            | 3.7 | 0.6 |
| TMS-XL7    | X-Large | 6.0                            | 3.7 | 0.7 |
| TMS-SIZ-NS | Sizers  | Set of 0.4cm Projection Sizers |     |     |

*Dimensions in centimeters*



### Conform™ Midfacial

With the added benefit of the soft and flexible Conform™ grid pattern, this midfacial implant augments both the malar and submalar region for a fuller, younger looking face. Sold in pairs.

| Cat. No.    | Size    | A            | B   | C    | D    |
|-------------|---------|--------------|-----|------|------|
| CCSM-XS     | X-Small | 3.4          | 2.0 | 0.39 | 0.20 |
| CCSM-S      | Small   | 4.1          | 2.5 | 0.47 | 0.27 |
| CCSM-M      | Medium  | 4.7          | 2.7 | 0.53 | 0.30 |
| CCSM-L      | Large   | 5.1          | 3.0 | 0.56 | 0.32 |
| CCSM-XL     | X-Large | 5.5          | 3.2 | 0.60 | 0.40 |
| CCSM-SIZ-NS | Sizers  | Complete Set |     |      |      |


*Dimensions in centimeters*

## Combined Submalar Shell™

Augments entire midfacial area by incorporating the features of the Terino Malar Shell® and the Binder Submalar®. Ideal for patients with both malar hypoplasia and midfacial soft-tissue deficiencies. Sold in pairs.

| Cat. No.   | Size    | A            | B    | C    | D    |
|------------|---------|--------------|------|------|------|
| CSM-S      | Small   | 4.38         | 3.00 | 0.40 | 0.30 |
| CSM-M      | Medium  | 4.90         | 3.26 | 0.50 | 0.40 |
| CSM-L      | Large   | 5.43         | 3.72 | 0.60 | 0.50 |
| CSM-XL     | X-Large | 5.78         | 4.16 | 0.70 | 0.60 |
| CSM-SIZ-NS | Sizers  | Complete Set |      |      |      |

*Dimensions in centimeters*




## Flowers Tear Trough®

Designed to augment the prominent suborbital groove caused by bony hypoplasia and depressions. A specified area of the implant is indented for removal to avoid the usual location of the infraorbital nerve. However, the actual position of the nerve should be verified to ensure the appropriate portion of the implant is excised. Sold in pairs.

| Cat. No.   | Size   | A            | B    | C    |
|------------|--------|--------------|------|------|
| FTT-S      | Small  | 4.22         | 2.32 | 0.25 |
| FTT-M      | Medium | 4.22         | 2.32 | 0.35 |
| FTT-L      | Large  | 4.75         | 2.60 | 0.45 |
| FTT-SIZ-NS | Sizers | Complete Set |      |      |

*Dimensions in centimeters*




## Extended Flowers Tear Trough®


Augments the prominent suborbital groove caused by bony hypoplasia and depressions, with continuous augmentation into the cheek area. Sold in pairs.

| Cat. No. | Size     | A   | B   | C    |
|----------|----------|-----|-----|------|
| EFTT-L   | Extended | 5.8 | 2.6 | 0.51 |

*Dimensions in centimeters*



## TEMPORAL IMPLANTS




### Temporal Shell

Designed by Barry L. Eppley, MD, DMD

An ideal solution for temporal hollowing due to genetics or trauma. The silicone Temporal Shell provides soft-tissue volume replacement to reestablish youthful convexity between the zygomatic process on up toward the forehead. Specify Catalog No. for "Left (L)" or "Right (R)". Single-packed (not in pairs). Save by purchasing a pair.

| Cat. No.  | A      | B            | C   |
|-----------|--------|--------------|-----|
| TS-04L    | 4.7    | 3.9          | 0.4 |
| TS-04R    | 4.7    | 3.9          | 0.4 |
| TS-07L    | 4.7    | 3.9          | 0.7 |
| TS-07R    | 4.7    | 3.9          | 0.7 |
| TS-09L    | 5.5    | 4.5          | 0.9 |
| TS-09R    | 5.5    | 4.5          | 0.9 |
| TS-SIZ-NS | Sizers | Complete Set |     |

*Dimensions in centimeters*



### Temporal Shell – Extended

Designed by Barry L. Eppley, MD, DMD

The extended temporal implant provides augmentation to both upper and lower temporal regions. From the top of the zygomatic arch to the anterior temporal line of the forehead, it provides correction of temporal hollowing through augmentation of the full vertical length of the muscle. Its tapered design allows the implant to blend smoothly into the temporal region. Specify Catalog No. for "Left (L)" or "Right (R)." Single-packed (not in pairs). Save by purchasing a pair.

| Cat. No.   | A      | B            | C   |
|------------|--------|--------------|-----|
| TS-04EL    | 4.0    | 7.0          | 0.4 |
| TS-04ER    | 4.0    | 7.0          | 0.4 |
| TS-06EL    | 4.5    | 7.5          | 0.6 |
| TS-06ER    | 4.5    | 7.5          | 0.6 |
| TSE-SIZ-NS | Sizers | Complete Set |     |


*Dimensions in centimeters*

## Flowers Dorsal Nasal

Provides mid-dorsal augmentation for central nasal depressions.

| Cat. No.   | Size   | A            | B    | C    |
|------------|--------|--------------|------|------|
| FDI-1.5    | Small  | 4.3          | 1.00 | 0.15 |
| FDI-2.5    | Sm-Med | 4.4          | 1.05 | 0.25 |
| FDI-3.5    | Medium | 4.4          | 1.10 | 0.35 |
| FDI-4.5    | Large  | 4.4          | 1.20 | 0.45 |
| FDI-SIZ-NS | Sizers | Complete Set |      |      |

*Dimensions in centimeters*




## Rizzo Dorsal Nasal

Provides full dorsal augmentation.

| Cat. No.   | Size    | A            | B    | C    | D    | E   | F   | G   |
|------------|---------|--------------|------|------|------|-----|-----|-----|
| RDI-XS     | X-Small | 4.5          | 2.25 | 0.20 | 0.15 | 0.5 | 0.7 | 0.9 |
| RDI-S      | Small   | 4.5          | 2.25 | 0.30 | 0.25 | 0.5 | 0.7 | 0.9 |
| RDI-M      | Medium  | 4.5          | 2.25 | 0.40 | 0.35 | 0.5 | 0.7 | 0.9 |
| RDI-L      | Large   | 4.5          | 2.25 | 0.47 | 0.45 | 0.5 | 0.7 | 1.0 |
| RDI-XL     | X-Large | 4.5          | 2.25 | 0.50 | 0.52 | 0.5 | 0.7 | 1.0 |
| RDI-SIZ-NS | Sizers  | Complete Set |      |      |      |     |     |     |

*Dimensions in centimeters*





## Peri-Pyriform™

A premaxillary implant designed to augment the prealar area around the pyriform aperture. Adds anterior projection to the premaxilla and provides a platform for the nasal base.

| Cat. No. | Size   | A   | B   | C   | D    |
|----------|--------|-----|-----|-----|------|
| BPI-S    | Small  | 3.8 | 1.5 | 2.8 | 0.53 |
| BPI-M    | Medium | 3.6 | 1.8 | 2.4 | 0.57 |
| BPI-L    | Large  | 3.9 | 2.1 | 2.8 | 0.65 |

*Dimensions in centimeters*






### Voloshin Dorsal Columella

Augments and reconstructs the nasal dorsum and the nasal tip. Fabricated in a soft and pliable opaque silicone elastomer.

| Cat. No.    | A      | B            | C   | D   | E   | F    | G    | H   |
|-------------|--------|--------------|-----|-----|-----|------|------|-----|
| VSDC-001    | 5.0    | 3.0          | 1.1 | 0.7 | 0.4 | 0.22 | 0.40 | 0.4 |
| VSDC-002    | 5.4    | 3.0          | 1.0 | 0.7 | 0.4 | 0.20 | 0.40 | 0.4 |
| VSDC-003    | 4.9    | 3.2          | 1.0 | 0.7 | 0.4 | 0.43 | 0.70 | 0.6 |
| VSDC-004    | 5.2    | 3.2          | 1.2 | 0.7 | 0.4 | 0.50 | 0.75 | 0.6 |
| VSDC-010    | 5.2    | 3.3          | 1.0 | 0.7 | 0.4 | 0.35 | 0.55 | 0.6 |
| VSDC-SIZ-NS | Sizers | Complete Set |     |     |     |      |      |     |

*Dimensions in centimeters*




### Anatomical Nasal Implant

Designed by Nazih M. Haddad, MD

Adds fullness to the dorsum. Provides a contoured tip and opens up the columella labial angle. Thin, tapered wings drape over the dorsal sides for implant stability and a gradual transition between the implant and surrounding tissue.

| Cat. No.   | A      | B            | C    | D    | E    | F    | G    | H    |
|------------|--------|--------------|------|------|------|------|------|------|
| ANI-010    | 3.7    | 2.3          | 0.80 | 0.66 | 0.80 | 0.20 | 0.58 | 0.42 |
| ANI-020    | 4.0    | 2.5          | 0.90 | 0.72 | 0.90 | 0.21 | 0.65 | 0.45 |
| ANI-030    | 4.4    | 2.7          | 1.00 | 0.79 | 0.95 | 0.23 | 0.70 | 0.50 |
| ANI-031    | 4.4    | 2.7          | 1.00 | 0.79 | 0.95 | 0.42 | 0.70 | 0.50 |
| ANI-040    | 4.7    | 2.9          | 1.05 | 0.84 | 1.00 | 0.23 | 0.75 | 0.55 |
| ANI-041    | 4.7    | 2.9          | 1.05 | 0.84 | 1.00 | 0.45 | 0.75 | 0.55 |
| ANI-SIZ-NS | Sizers | Complete Set |      |      |      |      |      |      |

*Dimensions in centimeters*



### Dorsal Columella

Carvable to fit each patient. Available in firm (F) and soft (S) durometers.

| Cat. No. | Type | A   | B   | C   | D   | E   | F   | G   |
|----------|------|-----|-----|-----|-----|-----|-----|-----|
| DCI-03F  | Firm | 5.8 | 3.0 | 1.1 | 0.3 | 0.9 | 0.8 | 0.6 |
| DCI-03S  | Soft | 5.8 | 3.0 | 1.1 | 0.3 | 0.9 | 0.8 | 0.6 |

*Dimensions in centimeters*

# SILICONE CARVING BLOCKS




Implantech's silicone carving blocks are easily fashioned into patient-specific implants. Choose from multiple shapes with varying durometers:

- Rectangular blocks
- Contoured blocks

## SILICONE CARVING BLOCKS



Clear



### Rectangular Carving Block – Clear

Carvable into individual implant shapes. Available in three hardnesses (durometers). Supplied non-sterile (sterilization instructions included).

| Cat. No.  | A    | B   | C   | Color | Durometer |
|-----------|------|-----|-----|-------|-----------|
| BLK-10-NS | 12.0 | 7.0 | 1.5 | Clear | Soft      |
| BLK-20-NS | 12.0 | 7.0 | 1.5 | Clear | Medium    |
| BLK-40-NS | 12.0 | 7.0 | 1.5 | Clear | Firm      |

*Dimensions in centimeters*



Opaque



Fleshtone


### Rectangular Carving Block – Various Durometers and Colors

Carvable into individual implant shapes. Available in five hardnesses (durometers): Very Soft (VS), Soft (S), Medium (M), Firm (F), Very Firm (VF) and four colors: Clear (C), Opaque (O), Fleshtone (F)

Supplied non-sterile (sterilization instructions included). May take 2-4 weeks processing time.

| Cat. No.    | A    | B   | C   | Color         | Durometer   |
|-------------|------|-----|-----|---------------|-------------|
| BL-1.0-C-NS | 13.6 | 5.8 | 1.0 | C (Clear)     | VS, S, M, F |
| BL-1.0-O-NS | 13.6 | 5.8 | 1.0 | O (Opaque)    | S, M, F     |
| BL-1.3-C-NS | 14.8 | 7.6 | 1.3 | C (Clear)     | VS, S, M, F |
| BL-1.3-O-NS | 14.8 | 7.6 | 1.3 | O (Opaque)    | S, M, F     |
| BL-2.0-C-NS | 14.9 | 7.7 | 2.0 | C (Clear)     | VS, S, M, F |
| BL-2.0-O-NS | 14.9 | 7.7 | 2.0 | O (Opaque)    | S, M, F     |
| BL-2.0-F-NS | 14.9 | 7.7 | 2.0 | F (Fleshtone) | S, M, F     |

*Dimensions in centimeters*




### Contoured Carving Block – Style 1 (Smooth)

Carvable and pre-formed in an elongated ovate shape with an asymmetrically tapered profile.

| Cat. No. | A    | B   | C   | Color | Durometer  | Volume |
|----------|------|-----|-----|-------|------------|--------|
| CCB1-1   | 11.9 | 5.3 | 1.5 | Clear | Extra Soft | 43cc   |
| CCB1-2   | 12.6 | 6.1 | 1.6 | Clear | Extra Soft | 58cc   |
| CCB1-3   | 13.7 | 6.8 | 1.7 | Clear | Extra Soft | 73cc   |
| CCB1-4   | 14.0 | 7.8 | 1.9 | Clear | Extra Soft | 102cc  |
| CCB1-5   | 16.2 | 9.4 | 2.0 | Clear | Extra Soft | 149cc  |

Individual sizes available

Dimensions in centimeters




### Contoured Carving Block – Style 3 (Smooth)


Carvable and pre-formed in a canoe shape with a symmetrically tapered profile.

| Cat. No. | A    | B   | C   | Color | Durometer  | Volume |
|----------|------|-----|-----|-------|------------|--------|
| CCB3-1   | 11.5 | 4.7 | 2.1 | Clear | Extra Soft | 55cc   |
| CCB3-2   | 14.1 | 5.0 | 1.8 | Clear | Extra Soft | 58cc   |
| CCB3-3   | 17.2 | 6.1 | 2.1 | Clear | Extra Soft | 93cc   |
| CCB3-4   | 20.3 | 6.5 | 2.1 | Clear | Extra Soft | 130cc  |

Individual sizes available

Dimensions in centimeters





### Contoured Carving Block – Style 4 (Smooth)


Carvable and pre-formed in an egg shape (ovoid).

| Cat. No.                   | A    | B   | Color | Durometer  | Volume |
|----------------------------|------|-----|-------|------------|--------|
| CCB4-1                     | 2.60 | 2.0 | Clear | Extra Soft | 5cc    |
| CCB4-2                     | 3.25 | 2.5 | Clear | Extra Soft | 10cc   |
| CCB4-3                     | 3.90 | 3.0 | Clear | Extra Soft | 19cc   |
| CCB4-4                     | 4.55 | 3.5 | Clear | Extra Soft | 29cc   |
| CCB4-5                     | 5.00 | 4.0 | Clear | Extra Soft | 40cc   |
| Individual sizes available |      |     |       |            |        |

*Dimensions in centimeters*

---


# 3D ACCUSCAN PATIENT-SPECIFIC® IMPLANTS



While our implant lines are quite comprehensive, stock items are occasionally insufficient for a patient's needs. As such, you can depend on Implantech's patient-specific products for an ideal solution.

Offerings include:

- 3D Accuscan® (for both the face and body)



### 3D Accuscan Patient-Specific® Implants

With 3D Accuscan®, FDA cleared for the face and entire body, you have the ability to quickly and easily reconstruct most bony defects and soft tissue deficiencies. Additionally, we can produce a 3D model of your patient's anatomy as a consultation or surgical aid. Contact us for additional information.

# SILICONE SHEETING




Implantech is the ideal source for all your silicone sheeting.

Choose from a variety of products to meet your medical applications, including:

- Silicone sheeting
  - ◆ Non-Reinforced
  - ◆ Long Term or Short Term Implantable

## SILICONE SHEETING - NON-REINFORCED




### AlliedSil™ Silicone Sheeting – Non-Reinforced, Long Term Implantable

Implantable for 30+ days. Supplied sterile. Box of 10 sheets.

| Cat. No.  | A | B | Thickness |
|-----------|---|---|-----------|
| 23-700-05 | 3 | 2 | .005      |
| 23-700-07 | 3 | 2 | .007      |
| 23-700-10 | 3 | 2 | .010      |
| 23-700-15 | 3 | 2 | .015      |
| 23-700-20 | 3 | 2 | .020      |
| 23-700-30 | 3 | 2 | .030      |
| 23-700-40 | 3 | 2 | .040      |

*Dimensions in inches*




### AlliedSil™ Silicone Sheeting – Non-Reinforced, Long Term Implantable

Implantable for 30+ days. Supplied sterile. Single sheet.

| Cat. No. | A | B | Thickness |
|----------|---|---|-----------|
| 700-05   | 8 | 6 | .005      |
| 700-07   | 8 | 6 | .007      |
| 700-10   | 8 | 6 | .010      |
| 700-15   | 8 | 6 | .015      |
| 700-20   | 8 | 6 | .020      |
| 700-30   | 8 | 6 | .030      |
| 700-40   | 8 | 6 | .040      |
| 700-50   | 8 | 6 | .050      |
| 700-60   | 8 | 6 | .060      |

*Dimensions in inches*



### AlliedSil™ Silicone Sheeting – Non-Reinforced, Short Term Implantable

Implantable for less than 30 days. Supplied sterile. Single sheet.

| Cat. No. | A | B | Thickness |
|----------|---|---|-----------|
| 600-05   | 8 | 6 | .005      |
| 600-07   | 8 | 6 | .007      |
| 600-10   | 8 | 6 | .010      |
| 600-15   | 8 | 6 | .015      |
| 600-20   | 8 | 6 | .020      |
| 600-30   | 8 | 6 | .030      |
| 600-40   | 8 | 6 | .040      |
| 600-50   | 8 | 6 | .050      |
| 600-60   | 8 | 6 | .060      |

*Dimensions in inches*

# SCAR MANAGEMENT



## ***Cimeosil®***

Cimeosil® Scar and Laser Gel & Cimeosil® Gel Sheeting are silicone-based products which have been shown to improve the appearance of scars. Cimeosil® has been used to reduce the visibility of scars from accidents, heart surgery, cesarean sections, abdominoplasty, breast reconstruction and many other causes.

## ***Gelzone® Wrap***

Made of stretchable, fabric-covered silicone gel, the patented Gelzone® Wrap combines compression, musculoskeletal support and scar management. The wrap stays secure while providing uniform compression. Beneficial for helping to flatten and fade keloids and hypertrophic scars.

## CIMEOSIL® SCAR AND LASER GEL & CIMEOSIL® GEL SHEETING




### Cimeosil® Scar and Laser Gel

Cimeosil® Scar and Laser Gel is topically applied to help flatten and fade scars and better manage scarring and redness from facial laser surgery. Made only of the highest quality medical-grade silicone, Cimeosil® Scar and Laser Gel is safe for both adults and children when used as directed. May be worn over sunscreen and under makeup. Choose from 5 gram and 14 gram tubes.

| Cat. No. | Description                               |
|----------|-------------------------------------------|
| CSG-05   | Cimeosil® Scar and Laser Gel-5 gram tube  |
| CSG-14   | Cimeosil® Scar and Laser Gel-14 gram tube |


*Recommended for older, deeper, more persistent scars. Gel sheeting is self-adhering, reusable and washable. All gel sheeting can be cut down to fit specific scar areas. A drying mat is included with each sheet.*



### Cimeosil® Gel Sheeting – Long Strip

To reduce large straight scarring resulting from cesarean sections, appendectomies, cardiac surgery, plastic surgery of the thighs and waist and non-surgical trauma. Sold individually, save by purchasing 6 units or more.



| Cat. No. | Description                     |
|----------|---------------------------------|
| GS-LS    | 4cm x 30cm Strip (1.6" x 11.8") |



### Cimeosil® Gel Sheeting – 4" x 5"

Multi-purpose - can be cut down to fit specific scar areas. Sold individually, save by purchasing 6 units or more.

| Cat. No. | Description   |
|----------|---------------|
| GS-04-05 | 4" x 5" Sheet |



### Cimeosil® Gel Sheeting – Areola Circle


To reduce scarring resulting from areola reconstruction or breast augmentation surgery. Sold individually, save by purchasing 6 units or more.

| Cat. No. | Description         |
|----------|---------------------|
| GS-AC    | Areola Circle Sheet |

### Cimeosil® Gel Sheeting – Vertical Mastopexy ‘Lollipop’

To reduce scarring from areola reconstruction, breast reduction or reconstructive surgery. Sold individually, save by purchasing 6 units or more.


| Cat. No. | Description              |
|----------|--------------------------|
| GS-VM    | Vertical Mastopexy Sheet |



### Cimeosil® Gel Sheeting – Breast Reduction ‘Anchor’

To reduce scarring of typical incisions that result from breast reduction or reconstructive surgery. Sold individually, save by purchasing 6 units or more.

| Cat. No. | Description            |
|----------|------------------------|
| GS-BR    | Breast Reduction Sheet |



### Cimeosil® Drying Mat

Each Cimeosil® Gel Sheeting product comes with a convenient drying mat. With its aerated grid, this plastic mat provides an ideal drying surface for your just-washed gel sheeting (not sold separately).



# GELZONE® WRAP



### Gelzone® Wrap

Made of stretchable, fabric-covered silicone gel, the patented Gelzone® Wrap is the only choice that combines compression healing, musculoskeletal support and scar management. The wrap stays exactly where it's put while providing ideal uniform compression. Excellent for flattening and fading keloids and hypertrophic scars.

*Benefits include non-slip compression for healing (after wound closure), complete silicone gel contact on closed wounds for scar management in difficult areas, general joint for muscle support and cold-wrap therapy. Hook-and-loop fastener tabs (extra tabs provided). All Gelzone® Wraps can be cut to size using scissors.*

### Gelzone® Wrap – 6" x 48"

Color: White  
Common Uses:  
Abdominal Support & Scar Management (After Liposuction); Abdominal Support & Scar Management (After C-Section); Abdominal Support & Scar Management (After Tummy Tuck); Arm Support & Scar Management (After Brachioplasty); General Joint or Muscle Support; Cold Wrap Therapy; Compression for Healing; Hamstring or Hip Pointer; Rib or Thorax Wrap; Osteoarthritis and Scar Management.



| Cat. No. | Description   |
|----------|---------------|
| WRAP-648 | 6" x 48" Wrap |

Thank you for your interest in Implantech's superior aesthetic products. We're pleased to offer various ways in which to order:

### **International Customers:**

Orders must be placed through an Implantech distributor in your country. To find a distributor, please access our "Distributor Locator" page at [www.implantech.com/distributor](http://www.implantech.com/distributor)

Once established with a distributor, you can access our website to view product information and educational resources.

Go on line at [www.implantech.com](http://www.implantech.com)

- New customers: Access "Int'l Physician Login" to register, then login
- Existing accounts: Access "Int'l Physician Login"

Call or email our international customer service department:

805-339-9415

[international@implantech.com](mailto:international@implantech.com)

### **Product Returns:**

All returns and exchanges must be pre-authorized. Please contact your local distributor to request a 'Return Goods Authorization' (RGA) number.

Patient-specific products and opened implant packages are non-returnable, and freight charges are ineligible for credit or refunds.

For our complete return/exchange policy, please refer to [www.implantech.com/return-policy](http://www.implantech.com/return-policy)





6025 Nicolle Street • Ventura, CA 93003  
800.733.0833 | [implantech.com](http://implantech.com)

